

therapy delivered  
**just where**  
needed

# CLINIPORATOR®

TECHNOLOGY FOR CARE



## ELECTROSCLEROTHERAPY **VASCULAR MALFORMATIONS**

## INTRODUCTION

Vascular malformations are congenital anomalies of the vascular system which can represent a therapeutic challenge, particularly in patients with recurrent or therapy-resistant malformations.


Current management mainly consists of intralesional sclerotherapy techniques. However, its effectiveness has been limited and often requires several treatment sessions.

**Electrosclerotherapy is a cutting-edge treatment for venous and lymphatic malformations.**


## ELECTROSCLEROTHERAPY

Electrosclerotherapy consists of the application of short-term, high intensity electrical pulses that induce transient permeabilisation of cell membranes (reversible electroporation), thus facilitating the intracellular delivery of an otherwise poorly permeant sclerosing agent.


Reversible electroporation induces the temporary disruption of cytoskeletal structures of endothelial cells and significantly enhances the effect of the sclerosing agent. Furthermore, electrical pulses have an important immediate vascular effect (vascular lock) that reduces blood flow transiently and prevents bleeding.<sup>[9-10]</sup>



## CLINICAL CASES



**Fig 1-A:** 10-year-old male patient. Symptomatic venous malformation of the left thigh unsuccessfully treated with surgery and sclerotherapy<sup>[2]</sup>



**Fig 1-B:** Complete resolution of the malformation three months after one electrosclerotherapy session<sup>[2]</sup>

## HOW EFFECTIVE IS ELECTROSCLEROTHERAPY

A recent observational study<sup>[2]</sup> conducted in 17 patients with therapy-resistant symptomatic venous malformations showed remarkable clinical effectiveness of electrosclerotherapy, with 8 patients becoming asymptomatic and 9 patients showing improvement after treatment.

These favourable results were achieved after only one electrosclerotherapy treatment session in 88% of the patients.

The median applied dose of bleomycin was 3 mg, which represents an additional advantage: lowering the risk of side effects related to a high cumulative dose.

The median volume reduction following electrosclerotherapy measured on MRI images was 86%.

**Electrosclerotherapy is used to treat different types of vascular malformations, including the following:**

- **Venous Malformations**
- **Macrocytic and microcystic lymphatic malformations**
- **Some forms of arteriovenous malformations**
- **FAVAs**

## BEYOND THE LIMITS WITH ELECTROSCLEROTHERAPY



**HIGHER  
EFFECTIVENESS**



**FEWER  
TREATMENT  
SESSIONS**



**REDUCED DOSE  
OF SCLEROSING  
AGENT**



**MINIMAL SIDE  
EFFECTS**



**Fig 2-A:** Therapy-resistant venous malformation



**Fig 2-B:** Result after 2 electrosclerotherapy sessions

Courtesy of Prof. W. Wohlgemuth  
University Clinic Halle (Germany)

# CLINIPORATOR®

TECHNOLOGY FOR CARE

## REFERENCES

1. Kostusiak M et al. Bleomycin Electrosclerotherapy Treatment in the Management of Vascular Malformations. *Dermatol Surg.* 2022 Jan 1; 48(1): 67-71.
2. Wohlgemuth WA et al. Bleomycin electrosclerotherapy in therapy-resistant venous malformations of the body. *J Vasc Surg Venous Lymphat Disord.* 2021 May; 9(3): 731-739.
3. Dalmády S et al. New Treatment Option for Capillary Lymphangioma: Bleomycin-Based Electrochemotherapy of an Infant. *Pediatrics.* 2020 Dec; 146(6): e20200566.
4. Horbach SER et al. Electrosclerotherapy as a novel treatment option for hypertrophic capillary malformations. *Dermatol Surg* 2020; 46: 491-8.
5. Gehl J et al. Updated standard operating procedure for cutaneous tumours and skin metastases. *Acta Oncol.* 2018 Jul;57(7):874-882.
6. Markelc B et al. Increased permeability of blood vessels after reversible electroporation is facilitated by alterations in endothelial cell-to-cell junctions. *J Control Release.* 2018 Apr 28; 276:30-41.
7. Mc Morrow L et al. Bleomycin electrosclerotherapy: new treatment to manage vascular malformations. *Br J Oral Maxillo fac Surg* 2017; 55: 977-9.
8. Horbach SER et al. Intralesional Bleomycin Injections for Vascular Malformations: A Systematic Review and Meta-Analysis. *Plast Reconstr Surg.* 2016 Jan; 137(1): 244-256.
9. Jarm T et al. Antivascular effects of electrochemotherapy: implications in treatment of bleeding metastases. *Expert Rev Anticancer Ther.* 2010 May; 10(5):729-46.
10. Kanthou C et al. The endothelial cytoskeleton as a target of electroporation-based therapies. *Mol Cancer Ther.* 2006 Dec; 5(12):3145-52.
11. Muir T et al. Bleomycin electrosclerotherapy (BEST) for the treatment of vascular malformations. An International Network for Sharing Practices on Electrochemotherapy (InspECT) study group report. *Radiol Oncol.* 2023 Jun 21;57(2): 141-149.
12. Guntau M et al. Individualized treatment of congenital vascular malformations of the tongue. *Clin Hemorheol Microcirc.* 2023; 83(4): 421-429.

### IGE A S.p.A.

HEADQUARTERS  
Via Parmenide, 10/A - 41012 Carpi (MO), ITALY  
Phone +39 059 699600  
[electroporation@igeamedical.com](mailto:electroporation@igeamedical.com)

### IGE A GERMANY

Feilitzschstr. 1,  
80802 München, GERMANY  
Phone +49 89 2304 1765  
[info.de@igeamedical.com](mailto:info.de@igeamedical.com)

### IGE A UK

Suites 1 & 2, Parkhill Business Centre Annexe,  
Walton Road, Wetherby, West Yorkshire LS22 5DZ, UK  
Phone: +44 1937 547065  
[info.uk@igeamedical.com](mailto:info.uk@igeamedical.com)



[www.igeamedical.com](http://www.igeamedical.com)